

# Increase Flexibility

For your daily liquid handling tasks with epMotion® 5070



## Applications

- > PCR/qPCR set-up
- > Normalization
- > Hit/Cherry picking
- > Dilution series
- > Plate reformatting

## Your benefits:

- > **More Space:** Higher capacity for plates and tips on deck on a small footprint
- > **More security:** Reduced pipetting errors from distractions
- > **Easy to use and flexible solution:** Use features like csv import on a large volume range to quickly start methods in your lab
- > **Gain time:** Automatic processing of up to 96 samples without user intervention
- > **Less risk:** Reduce the risk of RSI\* due to long pipetting series

\*repetitive strain injury

Get your demo and experience how to simplify your liquid handling

[www.eppendorf.com/epMotion5070](http://www.eppendorf.com/epMotion5070)

## What is new in epMotion 5070?

- > Modular rack to accommodate tips, reagents and/or consumables to free up deck capacity for e.g., more plates
- > Up to 3 deck positions for plates
- > Run full protocols without user-intervention
- > Fully automated processing of up to 96 samples for normalization, PCR or qPCR setup



## Ordering information

Description	Order no. International
<b>epMotion® 5070, EasyCon</b> , completely contained housing, system incl. EasyCon, epBlue™ software, mouse, waste box, 100 – 240 V ± 10% / 50 – 60 Hz ± 5%, 0.2 µL – 1 mL	5070006032
<b>epMotion® 5070 MultiCon</b> , completely contained housing, system incl. MultiCon, epBlue™ software, keyboard, mouse, waste box, 100 – 240 V ± 10% / 50 – 60 Hz ± 5%, 0.2 µL – 1 mL	5070000282
<b>Reservoir Rack Module PCR</b> , for PCR vessels and reagents	5075751933
<b>Reservoir Rack Module Tips</b> (set of 7 pieces, each tip module holds 16 tips)	5075751950

[www.eppendorf.com](http://www.eppendorf.com)

Eppendorf®, the Eppendorf Brand Design and epMotion® are registered trademarks of Eppendorf AG, Germany. epBlue™ is a trademark of Eppendorf AG, Germany. All rights reserved, including graphics and images. Copyright © 2019 by Eppendorf AG, Germany.