

# Statement on Safety

## Safety of Eppendorf mobile tables during operation with Eppendorf multipurpose centrifuges

The Eppendorf mobile tables

- Mobile table, low version (REF 5921 000.010)
- Mobile table, high version (REF 5921 000.012)

allow safe operation of the multipurpose centrifuges Centrifuge 5804, Centrifuge 5804 R, Centrifuge 5810, Centrifuge 5810 R, Centrifuge 5910 R, Centrifuge 5910 Ri, and Centrifuge 5920 R in compliance with the requirements of chapter 7.4 of the safety standard for laboratory equipment and laboratory centrifuges IEC 61010-1:2011 relating to IEC 61010-2-020:2007-03.

For the system consisting of Eppendorf mobile tables and Eppendorf multipurpose centrifuges of the above mentioned types, the same safety level is achieved as for operation on stationary laboratory workbenches.

The Eppendorf mobile tables are intended for use with the above specified Eppendorf multipurpose centrifuges.

The instructions for use of Eppendorf mobile tables must be read.

Date: March 2021



André Müller  
Product Life Cycle Manager Centrifugation



Monika Schneider  
Vice President Global Quality Management &  
Regulatory Affairs

**Your local distributor:** [www.eppendorf.com/contact](http://www.eppendorf.com/contact)  
Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany  
E-mail: [eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com)

ISO 9001  
Certified

ISO 13485  
Certified

ISO 14001  
Certified

Eppendorf® and the Eppendorf Brand Design are registered trademarks of Eppendorf AG, Germany. All rights reserved, incl. graphics and images. Copyright ©2021 by Eppendorf AG.